



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

## As per BIS 14543 Standard Approval Complete Lab List

### ANNEXURE 1 LIST OF EQUIPMENTS

| Sr. No. | Test equipment's/Chemicals And Identification numbers (Where Applicable)   | Least Count and Range (Where Applicable)   | Remarks (Indicate number of Equipment's) |
|---------|--|--|--|
| 1       | <b>Autoclave</b><br>Vertical Complete with pressure gauge, Safety Valve, Steam Release Valve Size 30 X 30 cm<br>Make :lab tech   | Range: 0 to 4.0. kg/ cm <sup>2</sup><br>Or<br>0 to 30 PSI<br><b>L.C- 0.05 kg/cm<sup>2</sup></b><br>Or<br>1 PSI | One                                      |
| 2       | <b>Hot air oven</b><br>size 12" X 12" X 12" made of MS with Thermostat control<br>Make :lab tech   | Temp 0 to 250 ° C<br>L.C: 2.0 ° C  | One                                      |
| 3       | <b>BOD Incubator 20° C</b><br>(With Digital Temp. Controller)<br>Size 14X14X14<br>Make : lab tech  | Temp 0 to 50 ° C<br>L.C: 0.1 ° C   | one                                      |
| 4       | <b>BOD Incubator 25° C</b><br>(With Digital Temp. Controller)<br>Size 12"<br>Make : lab tech   | Temp 0 to 50 ° C<br>L.C: 0.1 ° C   | one                                      |
| 5       | <b>Incubator Digital temp 44°C</b><br>bacteriological<br>size 14" X 14" X 14" inner chamber of SS outer of MS with thermostat (Thermostatically controlled With Thermometer)<br>Make :lab tech |  | one                                      |
| 6       | <b>Incubator Digita temp 37°C</b><br>bacteriological<br>size 14" X 14" X 14" inner chamber of SS outer of MS with thermostat (Thermostatically controlled With Thermometer)<br>Make : lab tech |  | one                                      |
| 7       | <b>Gas Burner</b>  | -  | One                                      |
| 8       | <b>Water bath (With glass Thermometer)</b><br>inner chamber of SS, outer of SS 6 holes of 3" with Digital control complete SS<br>Make : lab tech   | Range: 0 to 110 ° C<br>LC : 1.0 ° C  | One                                      |
| 9       | <b>Laminar air flow</b>  | 1.5 feet UV Tube   | One                                      |
| 10      | <b>Hot plate</b>   | Range: 0 to 100 ° C<br>L.C- 1 ° C  | One                                      |



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

| Sr. No. | Test equipment's/Chemicals And Identification numbers (Where Applicable)              | Least Count and Range (Where Applicable)                | Remarks (Indicate number of |
|---------|---|---|-----------------------------|
| 11      | <b>pH meter digital</b><br>Make :INSIF INDIA  | 0-14 pH, LC- 0.01 pH                                    | One                         |
| 12      | <b>spectrophotometer</b><br>Digital resolution 1 nm<br>Make : INSIF INDIA             | 360 to 960 nm<br>LC- 1 nm % t,<br>0.001 ABS, 0.001 Cons | One                         |
| 13      | <b>Turbidity counter digital</b><br>Make : INSIF INDIA                                | 0 to 1000 NTU<br>LC : 0.1 NTU                           | One                         |
| 14      | <b>Colony counter digital</b><br>Make : chemi line                                    | Range :0 to 999 Count<br>L.C- 1 Count                   | One                         |
| 15      | <b>S.S. Filter holder assembly</b>  | 47 mm   | One                         |
| 16      | <b>Microscope Binocular</b>   | 10 X, 15X,  | One                         |
| 17      | <b>Magnetic stirrer with hot plate</b>  | Cap : 1 Litre<br>Range : 0 to 100%<br>L.C- 1 %          | One                         |
| 18      | <b>Centrifuge</b>   | 4 X15 ml  | One                         |
| 19      | <b>Single pan balance</b><br>Make :U-MAX  | <b>0 to 200 gm</b><br><b>L.C : 0.1 gm</b>               | One                         |
| 20      | <b>High pressure vacuum pump :</b>  | ¼ HP motor<br>50 lit per minute                         | One                         |
| 21      | <b>Gas stripping apparatus</b><br>Borosilicate glass with fuming hood & Air flow pump | 1 litre   | One                         |
| 22      | <b>Flow meter</b>   | Range: 0 to 100 LPH<br>L.C : 5 LPH                      |                             |
| 23      | <b>Heating mentalis</b><br>Make :---  | 2 litre<br>range: 0 to 100%<br>L.C- 1 %                 | One                         |
| 24      | <b>Florescent Inspection Screen</b>   | UV Light 9" long  | One                         |
| 24      | <b>Carbon Column</b>  | 150 mm<br>250 ml  | one                         |
| 25      | <b>Distillation Assembly</b>  | 2 Ltr / Hr  | one                         |



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : +91 73595 95577  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : +91 70432 55777  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : +91 281 - 2922777  
Telefax : +91 281 - 2922481  
Cell : +91 74050 50444 / +91 70439 55777  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

### ANNEXURE 2 LIST OF Chemical

| Sr. No. | Chemicals and its grade<br>(Where Applicable) | Pack    |
|---------|---|---------|
| 01      | Activated charcoal                            | 100gm   |
| 02      | Acetic acid                                   | 500ml   |
| 03      | Aluminum oxide                                | 100gm   |
| 04      | Aluminum potassium sulphate Ar                | 100gm   |
| 05      | Ammonium Acetate                              | 100gm   |
| 06      | Ammonium Carbonate                            | 100gm   |
| 07      | Ammonium Di chromate                          | 100gm   |
| 08      | Ammonium Ferrous sulphate                     | 100gm   |
| 09      | Ammonium metavanadate                         | 5gm     |
| 10      | Ammonium persulphate AR                       | 5gm     |
| 11      | Di Ammonium Hydrogen orthophosphate           | 100gm   |
| 12      | Ammonium Thiocyanate                          | 100gm   |
| 13      | Ammonia solution                              | 500ml   |
| 14      | Antimony metal powder                         | 25gm    |
| 15      | Ascorbic acid                                 | 25gm    |
| 16      | Barium chloride                               | 100gm   |
| 17      | Barium Hydroxide                              | 100gm   |
| 18      | Bromo cresol green                            | 100ml   |
| 19      | Boric Acid                                    | 100gm   |
| 20      | Buffer tablet P.H.4.0                         | P.H.4.0 |
| 21      | Buffer tablet P.H.7.0                         | P.H.7.0 |
| 22      | Calcium Carbonate                             | 100gm   |
| 23      | Calcium Hydroxide                             | 100gm   |
| 24      | Chromotropic acid                             | 5gm     |
| 25      | Citric Acid                                   | 100gm   |
| 26      | Cupric sulphate                               | 100gm   |
| 27      | Chloroform                                    | 500ml   |
| 28      | Cobalt chloride                               | 25gm    |
| 29      | Copper metal powder                           | 25gm    |
| 30      | Cyclohexanone                                 | 500ml   |
| 31      | EDTA Di Sodium salt                           | 100gm   |
| 32      | Di ethylene glycol monobutyl ether            | 500ml   |
| 33      | Di (2 ethyl hexyl) sulphosuccinate            | 25gm    |
| 34      | Dodecyl benzene sulphonic acid                | 100gm   |
| 35      | Eriochrome Black T                            | 25gm    |
| 36      | Eriochrome cyanine R AR                       | 5gm     |
| 37      | Ethylene glycol                               | 500ml   |
| 38      | Ethyl acetate                                 | 500ml   |
| 39      | Ferric Chloride                               | 100gm   |
| 40      | Ferric Citrate                                | 100gm   |



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

| Sr. No. | Chemicals and its grade<br>(Where Applicable) | Pack  |
|---------|---|-------|
| 41      | Formaldehyde                                  | 500ml |
| 42      | Gelatin                                       | 100gm |
| 43      | Glycerine                                     | 500ml |
| 44      | Hexamine                                      | 100gm |
| 45      | Hydrochloric acid                             | 500ml |
| 46      | Hydrogen peroxide 30%                         | 250ml |
| 47      | Hydroxylamine hydrochloride                   | 25g m |
| 48      | Di iso propyl ether                           | 100ml |
| 49      | Hydrazine sulphate                            | 100gm |
| 50      | Iso propyl alcohol                            | 500ml |
| 51      | Manganase Sulphate                            | 100gm |
| 52      | Methanol                                      | 500ml |
| 53      | Mercuric chloride                             | 2 gm  |
| 54      | Methylene blue                                | 100ml |
| 55      | Methyl Orange                                 | 100ml |
| 56      | Methyl Red                                    | 100ml |
| 57      | Nitric acid LR                                | 500ml |
| 58      | NED Di hydrochloride                          | 5 gm  |
| 59      | Neocuprine AR                                 | 1gm   |
| 60      | N,N,Di methyl P-Phenylene Di amine oxalate AR | 2gm   |
| 61      | Orthophosphoric acid                          | 250ml |
| 62      | Patton & Readers                              | 5 gm  |
| 63      | 1,10 phenanthroline                           | 5 gm  |
| 64      | Phenophtheline                                | 100ml |
| 65      | Potassium Hydrogen sulphate                   | 100gm |
| 66      | Potassium chloride AR                         | 100gm |
| 67      | Potassium Chromate                            | 100gm |
| 68      | Potassium Dichromate                          | 100gm |
| 69      | Potassium Iodide                              | 25gm  |
| 70      | Potassium ferricyanide                        | 100gm |
| 71      | Potassium chloroplatinate                     | 100ml |
| 72      | Potassium di hydrogen phosphate               | 100ml |
| 73      | Di potassium hydrogen phosphate               | 100gm |
| 74      | Potassium hydroxide                           | 100gm |
| 75      | Potassium meta periodate                      | 25gm  |
| 76      | Potassium nitrate                             | 100gm |
| 77      | Potassium permanganate                        | 100gm |
| 78      | Silica gel blue                               | 100gm |
| 79      | Silver nitrate                                | 5 gm  |
| 80      | Sodium Acetate                                | 100gm |



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

| Sr. No. | Chemicals and its grade<br>(Where Applicable) | Pack  |
|---------|---|-------|
| 81      | Sodium arsenite                               | 5gm   |
| 82      | Sodium ascorbate                              | 25gm  |
| 83      | Sodium bi carbonate                           | 100gm |
| 84      | Sodium carbonate                              | 100gm |
| 85      | Sodium chloride                               | 100gm |
| 86      | Sodium citrate                                | 100gm |
| 87      | Sodium hydroxide                              | 100gm |
| 88      | Sodium hypochloride                           | 500ml |
| 89      | Sodium nitrate                                | 100gm |
| 90      | Sodium nitrite                                | 100gm |
| 91      | Sodium oxalate                                | 100gm |
| 92      | Sodium sulphate                               | 100gm |
| 93      | Sodium sulphite                               | 100gm |
| 94      | Sodium thiosulphate                           | 100gm |
| 95      | Sulphur dioxide solution                      | 100ml |
| 96      | Sulphuric acid AR                             | 500ml |
| 97      | Sulphanilamide                                | 25gm  |
| 98      | Starch  | 100gm |
| 99      | Triethanolamine                               | 500ml |
| 100     | Urea  | 100gm |
| 101     | Zinc acetate                                  | 100gm |
| 102     | Zinc sulphate                                 | 100gm |
| 103     | Zinc amalgamated                              | 10gm  |
| 104     | Zinc dust AR                                  | 25gm  |
| 105     | Zincon AR                                     | 1gm   |





MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

### Microbiology Media

|     |  |        |
|-----|--|--------|
| 1.  | Plate count agar   | 100gm  |
| 2.  | Asparagine (DL) proline (L) broth  | 100gm  |
| 3.  | Differential clostridium reinforce broth   | 100gm  |
| 4.  | EMB agar   | 100gm  |
| 5.  | chloramphenicol Subouroud Agar   | 100gm  |
| 6.  | MacKonkey's agar   | 100 gm |
| 7.  | Lactose TTC agar with sodium heptadecylsulphat (lactose, peotone, yeast extract, meat extract, bromothymol blue, agar: 2,3,5, triphenyltrazolium chloride (TTC), sodium heptadecylsulphat (tergistoal-7) | 100 gm |
| 8.  | Tryptophan broth (tryptic digest of casine, L-tryptophan, sodium chloride)   | 100gm  |
| 9.  | Tryptone soy agar (TSA)-(tryptic digest of casine, soy peptone, NaCL, agar)  | 100 gm |
| 10. | Tryptone bile agar (TBA)   | 100 gm |
| 11. | Oxidase reagent (tetramethyl-p-phenylene diamine hydrochloride)  | 1 GM   |



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

## ANNEXURE 3

GLASS WARE-1

| Sr. No. | Description               | Liter / ML | Qty.    |
|---------|---------------------------|------------|---------|
| 1       | Nessler Cylinder          | 50 ML      | 5       |
| 2       | Standard Volumetric Flask | 10 ML      | 4       |
|         |                           | 50 ML      | 4       |
|         |                           | 100 ML     | 4       |
|         |                           | 250 ML     | 2       |
|         |                           | 500 ML     | 2       |
|         |                           | 1000 ML    | 1       |
| 3       | Graduated Pipette         | 1 ML       | 1       |
|         |                           | 2 ML       | 1       |
|         |                           | 5 ML       | 1       |
|         |                           | 10 ML      | 1       |
|         |                           | 25 ML      | 1       |
| 4       | Separating Funnel         | 125 ML     | 1       |
|         |                           | 500 ML     | 1       |
| 5       | Conical Flask             | 250 ML     | 8       |
|         |                           | 500 ML     | 2       |
| 6       | Beaker (Glass)            | 50 ML      | 3       |
|         |                           | 100 ML     | 3       |
|         |                           | 250 ML     | 3       |
|         |                           | 500 ML     | 3       |
| 7       | Petri Dish                | -----      | 20 Pair |
| 8       | Screw Cap Bottle          | 30 ML      | 4       |
|         |                           | 100 ML     | 4       |
|         |                           | 250 ML     | 2       |
| 9       | Iodine Flask              | 100 ML     | 1       |
| 10      | Filtration Flask          | 1000 ML    | 1       |
| 11      | Carbon Column             | 1 Feet     | 1       |
| 12      | Wide Mouth Bottle         | 1000 ML    | 1       |



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

## ANNEXURE 4

GLASS WARE-2


| Sr. No. | Description                | Liter / ML | Qty. |
|---------|----------------------------|------------|------|
| 1       | Beaker (Plastic)           | 50 ML      | 2    |
|         |                            | 100 ML     | 2    |
|         |                            | 250 ML     | 2    |
|         |                            | 500 ML     | 2    |
| 2       | Wash Bottle                | 250 ML     | 2    |
| 3       | Burette Stand              | -----      | 1    |
| 4       | Separating Funnel Stand    | -----      | 1    |
| 5       | Funnel                     | -----      | 2    |
| 6       | Stand for Nessler Cylinder | -----      | 1    |
| 7       | Pipette Stand              | -----      | 1    |
| 8       | Measuring Cylinder         | 50 ML      | 1    |
|         |                            | 100 ML     | 2    |
| 9       | Spatula                    | 150 MM     | 1    |
| 10      | Reagent Bottle (Plastic)   | 250 ML     | 24   |
|         |                            | 500 ML     | 12   |
| 11      | Glass Slide                |            | 1    |
| 12      | Burner                     |            | 2    |
| 13      | Inoculation loop           |            | 1    |
| 14      | Pipette Bulb               | Small      | 1    |
|         |                            | Medium     | 1    |
|         |                            | Large      | 1    |
| 15      | Tissue Paper               |            | 1    |
| 16      | Cotton                     |            | 1    |
| 17      | Desiccator                 |            | 1    |



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)



### Additional Arrangements to be done By Client:

1. Ac. in Microbiology Lab.
2. LPG gas Cylinder with Regulator and tuning for Micro Lab.
3. Washing Powder and Liquid soap in both labs.
4. Ventilation fan in Chem. lab.
5. Gas Stripping app chamber- 1.5 'depth, 2' width, 3 ' height
6. Electric plugs: 10 in Chem. lab and 11 in Microbial Lab.
7. Height of platform : 3feet, 2.5 to 2.8 feet out
8. Geezer, Mask, Hand Gloves & Apron of Workers
9. First Aid Box for 10 Persons Insect Killer for two sites.
10. 14543 ISI testing standard.
11. Refrigerator Videocon 100 liter (1-nos.).

### Extra Additional to be done By Client Amount Rs/-:

- 123000/- -BIS License Fees Yearly.
- 12000/- -Row water taste report.
- 9440/- -BIS application Fees.
- 10000/- -BARC report.
- 50000/- -Sample testing fees(jar, bottle, cup).
- 12000/- -ISI standard 14543 book.
- 25000/- to 50000/- -BIS officer.
- 25000/- For BIS File Work



MUMBAI OFFICE :  
116-D, Govt. Industrial Estate,  
Charkop, Kandiwali (West),  
Mumbai - 400067,  
Maharashtra.  
Cell : [+91 73595 95577](tel:+917359595577)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

AHMEDABAD OFFICE :  
B/9, Sumel-7, Business Park,  
Opp. Soni ni Chal,  
BRTS Bus Stop, N.H.-8  
Ahmedabad - 382325.  
Cell : [+91 70432 55777](tel:+917043255777)  
E-mail : [inquiry@rajwater.com](mailto:inquiry@rajwater.com)  
[info@rajwater.com](mailto:info@rajwater.com)

HEAD OFFICE :  
Plot No. : 27, Survey No. 47, Jivraj Industrial Area,  
Nr. Falcon Pump, Vavdi, Gondal Road, Rajkot - 360004.  
Tel : [+91 281 - 2922777](tel:+912812922777)  
Telefax : [+91 281 - 2922481](tel:+912812922481)  
Cell : [+91 74050 50444](tel:+917405050444) / [+91 70439 55777](tel:+917043955777)  
Website : [www.rajwater.com](http://www.rajwater.com)  
Email : [inquiry@rajwater.com](mailto:inquiry@rajwater.com) / [info@rajwater.com](mailto:info@rajwater.com)

